

Dry Biochemical Analyzer



Low-maintenance



High precision



Disposable consumable



≤3min



4.3 Inch
touch screen



Constant temperature
and humidity



Ultra portable



Optional printer

Test Parameters:

Lipids + Glucose

TC | TG | HDL-C | LDL-C | GLU

Kidney Function

BUN | CRE | UA

Liver Function

ALB | ALT | AST | ALP


Screening for Metabolic diseases

TC | UA | GLU


Screening of blood donors

Hb | ALT


Testing in three easy steps



1.puncture the finger using a lancet



2.draw a sufficient amount of blood when the blood drop is large enough



3.drop the blood into the sample area of the test card for testing

Specifications:

Detection method	Photometric detection	Measurement	157.5*145.7*64.6mm
Detection time	≤ 3min	Weight	622g
CV	≤ ± 10%	Data storage	3000 results
Sample type	Peripheral or venous blood or serum	Display	4.3touch screen
Sample volume	≤ 45 μL	Storage environment	Temp:20℃ ~60℃ ; Rh:10% ~ 90%
Power supply	Adapter input AC100~240V,50/60Hz,30VA; output5V,2A Built-in Lithium battery: DC3.7±0.3V	Operation environment	Temp:15℃ ~35℃ ; Rh:≤85%

JIANGSU KONSUNG BIO-MEDICAL SCIENCE AND TECHNOLOGY CO.,LTD

Email: info@konsung.com

Global marketing center: Room 606, Building 2, E-PARK, 47 Andermen Street, Nanjing, China

R & D Center: A01, 14th Floor, Building 1, Phase 2, New Strategic Emerging Industrial Park, Liuxian Road, Shenzhen city, China

Produce center: No.8, West Shengchang Rd, Danyang Economic Development zone, Danyang city, Jiangsu, China

Website: <http://en.konsung.com/>