

Product Manual

HEMATOCRIT CENTRIFUGE HC - 12C

- ❖ With brake and automatic alarm after operation.
- ❖ With mechanical door lock and safety switch.
- ❖ Strong Plastic case with metal lid.
- ❖ High power brushless motor
- ❖ Over speed and imbalance alarm
- ❖ With Reader.
- ❖ Accurately measuring the specific value of blood cell.
- ❖ Separating micro-amount of blood or solution.



MODEL	HC-12C
Max RPM	12000 RPM. Stepless speed adjustment, LED display speed and time
MAX RCF	15300xg
TIME	0-99 min
Rotor Capacity	Ø 1.5x75mmx24 capillary tubes
Voltage	110V 60Hz/ 220V 50Hz
Dimensions (mm)	280x310x265
Weight	10 kg

Zenith Lab Inc.
1979 W Holt Ave, Pomona,
CA 91768 USA
Tel: 1-909-868-1800
Fax: 1-909-868-1898
Email: info@zenithlabo.com

