

冷却小型卓上遠心機

H-19FMR



仕様 SPECIFICATIONS

	H-19FMR
Max. Speed	4,300rpm
Max. R.C.F	3,000×g
Max. Capacity	400ml
Control	Microcomputer
No. of Memories	3 Memories
Power-Source	AC100V 1φ 5A 1KVA
Dimensions	W430×D460×H400mm
Weight	40kg
Note	

小型卓上遠心機

H-19α



仕様 SPECIFICATIONS

	H-19α
Max. Speed	4,000rpm
Max. R.C.F	2,650×g
Max. Capacity	400ml
Control	Microcomputer+Inverter
No. of Memories	10 Memories
Power-Source	AC100V 1φ 5A 0.5KVA
Dimensions	W364×D384×H320mm
Weight	21.5kg
Note	

H-19Fシリーズ適用ロータ

Type	Capacity	Note
RF-109K	Rotor Swing	
MA-109K	15ml×4 Swing 4,000rpm 2,600×g	
MB-109K	15ml×8 Swing 4,000rpm 2,600×g	
MC-109K	15ml×16 Swing 4,000rpm 2,600×g	
MD-109K	50ml×4 Swing 4,000rpm 2,600×g	
ME-109K	50ml×4(TC) Swing 4,300rpm 3,000×g	✓
MF-109K	15ml×24 Swing 4,000rpm 2,600×g	✓
PA-109	10ml×24 Swing 3,000rpm 1,460×g	
PB-109	15ml×16 Swing 3,000rpm 1,430×g	
PC-109	6ml×40 Swing 3,000rpm 1,160×g	
PD-109	1.5(2.0)ml×48 Swing 3,000rpm 990×g	
PE-109	10ml×16 Swing 3,000rpm 1,460×g	
PN-109	15ml×24 Swing 3,500rpm 1,950×g	✓
PL-109	5ml採血管×24 Swing 3,000rpm 1,160×g	
RF-109L	Rotor Swing	
MF-109L	100ml×4 Swing 4,000rpm 2,650×g	
RA-011	15ml×12 Angle 4,000rpm 2,040×g	
RA-012	15ml×8 Angle 4,000rpm 2,040×g	
RA-013	10ml×16 Angle 4,000rpm 2,040×g	
RA-014	7ml×24 Angle 4,000rpm 2,040×g	
RA-015	50ml×4 Angle 3,500rpm 1,740×g	

H-19FMRは、(H-19F)H-19FMR(4,000rpm)



RF-109K



MA-109K



MB-109K



MC-109K



MD-109K



ME-109K



MF-109K



PF-109L



PA-109



PB-109



PC-109



PD-109



PE-109



PN-109



PL-109