

Urine Chemistry Analyzer


Reader 300

New

Ref No. : 5200

Product Features

- ◆ User friendly designed
- ◆ Easy to use and efficient working cycle with multi-strip tray
- ◆ Optional operation modes
(General / One-by-One / Quick mode)
- ◆ High speed throughput
- ◆ Compatible with various kinds of CYBOW strips
- ◆ Easy to clean and do maintenance
- ◆ Program update by S/W downloading
- ◆ Multi-languages are available
- ◆ Accessible to Laboratory Information System (LIS)



Technical Specification

Type : Reflectance Photometer	Throughput : 300 Tests / hour (Max 1,000 Tests/hour)
Dimension : 275(W) × 250(D) × 170(H) mm	Interface : RS232c
Net Weight : 1.3 kg	ID Input: Barcode Scanner or Keyboard
Power Supply : 12V DC / 3 A	Memory : 2,000 Tests
Printer : Thermal printer	Compatible strips : CYBOW 10M, 11M, 2MAC, 12MAC etc.
Wave length : 460, 550, 650nm	Temperature (°C) / Humidity (%) : 2~30 °C / 10~70 %

Basic Accessories: AC Adaptor, Power Cable, Operation Manual, RS232C cable, Calibration strips

Optional Accessories: External Keyboard, Barcode Scanner

Product Certificates



542-1, Daman-Ri, Jinrye-Myun, Gimhae-City
621-881, Gyung-Nam, Korea (South)
Tel: 82-55-3461882 / Fax: 82-55-3461883
E-mail: sales@cybow.com
Web-Site: www.cybow.com

(Rev. - 04)