



Urine Chemistry Analyzer

Reader 720

Ref No. : 5400

Product Features

- ◆ Semi-Automatic Analyzer
- ◆ designed for multiple tests
- ◆ Compact size with High Throughput
(Max 720 Tests per Hour)
- ◆ Automatic start of measurement
after placement of strip
- ◆ Easy to clean and do maintenance
- ◆ Program update by S/W downloading
- ◆ Multi-languages are available
- ◆ Accessible to Laboratory Information System (LIS)
- ◆ Compatible with various kinds of CYBOW strips



Technical Specification

Type: Reflectance Photometer	Throughput: 720 Tests / hour
Dimension: 320(W) x 260(D) x 178(H) mm	Interface: RS232C
Net Weight: 2.4 kg	ID Input: via Barcode Reader
Power Supply: 12V DC / 3A	Memory: 2,000 Tests
Printer: Thermal Printer	Compatible strips: CYBOW 10M, 11M, 2MAC, 12MAC etc.
Wave length : 460, 550, 650nm	Temperature (°C) / Humidity (%) : 2~30 °C / 10~70 %

Basic Accessories: AC Adaptor, Power Cable, Operation Manual, RS232C cable, Calibration strips

Optional Accessories: External Barcode Scanner

Product Certificates



542-1, Daman-Ri, Jinrye-Myun, Gimhae-City
621-881, Gyung-Nam, Korea (South)
Tel: 82-55-3461882 / Fax: 82-55-3461883
E-mail: sales@cybow.com
Web-Site: www.cybow.com

(Rev. - 04)