


CYBOW™ R-50

Urine Chemistry Analyzer



Ref No. : 5100

Product Features

- ◆ Semi-Automatic Analyzer Designed For One-by-One Test
- ◆ Built-in Printer and Wide LCD Monitor
- ◆ Easy and Convenient Operation with Shortcut Buttons
- ◆ Compatible with Various kinds of CYBOW Strips
- ◆ Accessible to Laboratory Information System (LIS) by PC protocol

Technical Specification

Type: Reflectance Photometer	Throughput: 45 Tests / hour (Max 120 Tests/hour)
Dimension: 300(W) x 250(D) x 90(H) mm	Interface: RS232C
Net Weight: 1.2 kg	ID Input: via External Keyboard
Power Supply: 12V DC / 2A	Memory: 2,000 Tests
Printer: Thermal Printer	Compatible strips: CYBOW 10, 11, 2AC, 10AC

Basic Accessories: AC Adaptor, Power Cable, Operation Manual, RS232C cable, Calibration strips

Optional Accessories: External Keyboard

■ Product Certificates



542-1, Daman-Ri, Jinrye-Myun, Gimhae-City
621-881, Gyung-Nam, Korea (South)
Tel: 82-55-3461882, 3452593 / Fax: 82-55-3461883
E-mail: sales@cybow.com
Web-Site: www.cybow.com